

Medit SmartX Integration Partners

Last updated: October 22, 2025

Integrated Partner- Product Name	Logo	Available Market*	System Compatibility**	Medit SmartX 2.0 Integration** (Anchor Cap Support)
Scan Ladder - Titanium Direct Kit		US  EU  UK  JP  AU  + Global	-	0 / Tissue Sync
Apollo Implant Components - SmartFlag		US  EU 		X
Dentis - SCANBODY PRO		US  EU  KR  IN 		X
GeoMedi - ArchB Scanbody Solution		KR  JP 		X
Segma - Digital Wings		CN  US  EU 		X
Top implant - Truefit Pro3		CN  US  EU 		X
Digital Arches - OptiSplint		US 		0 / SCANCAPS
Implant Direct - DIRECTscan		US 		X
RR Dental Lab - FireX		IN 		X
NeoBiotech - SmileArch ScanBody Kit		US  KR 		X
Derbarzt - Derbarzt Scan System		EU  KR 		0 / Gingiva Scanbody

* Available Market

The registration status and list of available markets may change over time. For the most up-to-date information, please refer to each manufacturer's official website.

** System Compatibility

All partner solutions are compatible with the industry-standard **Nobel system**, ensuring seamless interoperability. In addition, each manufacturer provides extended compatibility with various implant systems, offering users greater flexibility and choice. For detailed compatibility information, please refer to each manufacturer's official website or product documentation.

*** What is an Anchor Cap?

An Anchor Cap is a cap-shaped digital library automatically generated based on the implant position after scan body alignment. This feature allows clinicians to utilize both **horizontal (archform) scan body** workflows and **centralized scan body** workflows, expanding flexibility and compatibility across a wide range of clinical applications.

When the Anchor Cap is physically attached and scanned inside the patient's mouth, it serves as a stable reference point (anchor) for soft tissue scanning. This enhances overall scan accuracy and ensures smooth, reliable data capture even in challenging cases such as immediate loading, where bleeding or tissue deformation can affect recognition.

 Note: The Anchor Cap is a standardized concept defined within Medit SmartX 2.0, but actual product names may vary depending on the manufacturer. For example: Scan Ladder: *Tissue Sync*, Optisplint: *SCANCAPS*, Derbarzt: *Gingiva Scanbody*